

# Certificate of Compliance

This is to certify that the enclosed check source(s) conforms to the conditions and limitations specified in the regulations listed below.

- **U.S. NRC – 10 CFR Part 30.18 “Exempt Quantities”**

**NOTE:** The source(s) in this package does not exceed the activity limit set forth in Section 30.71 Schedule B of the U.S. NRC Regulations (the source may be transferred to any individual in the U.S. free of specific licensing requirements); however, it may be licensable material in countries other than the U.S.

- **U.S. DOT – 49 CFR Part 173.421 "Excepted Packages for Limited Quantities of Class 7 (Radioactive) Materials"**

- **USPS – Postal Publication No. 52 "Hazardous, Restricted, and Perishable Mail"**

- **IAEA – No. TS-R-1 "Regulations for the Safe Transport of Radioactive Materials"**

- **IATA – Section 4.2 "Identification" (UN2910 Radioactive Material, Excepted Package, Limited Limited Quantity of Material N.O.S., Class 7)**

- **IATA – Section 10.3.11.1.2. "Radioactive Material in Limited Quantities"**

**All sealed sources are wipe tested for surface contamination prior to shipping.**

**NOTE:** Unsealed sources (e.g. needle sources, Po-210 disk sources), which have an active area that is uncovered or protected by a very thin coat, are not designed to pass a standard leak test even though the deposit is adherent. The inactive portions of the source have been checked using a standard wipe test.

**Product Description:**

Model No.	Nuclide	Nominal Activity	Half-life	Quantity
-----------	---------	------------------	-----------	----------

Authorized by: \_\_\_\_\_

Andrew Cordell, Radioisotope Manager

Date: \_\_\_\_\_

**Manufactured by:**

**SPECTRUM TECHNIQUES, LLC.**

106 Union Valley Road, Oak Ridge, TN 37830

Phone: 865.482.9937 Fax: 865.483.0473 Email: sales@spectrumtechniques.com

www.spectrumtechniques.com